

## Application Guide for Viscometer

Consider application and viscosity range when selecting model (LV,RV,HA,HB)

### LV SERIES - LOW VISCOSITY

Adhesives (solvent base)	Juices	Photo Resist
Biological Fluids	Latex	Polymer Solutions
Chemicals	Milk	Rubber Solutions
Dairy Products	Oils	Solvents
Hot Waxes	Paints and Coatings	
Inks	Pharmaceuticals	

### RV SERIES - MEDIUM VISCOSITY

Adhesives (hot melt)	Food Products	Paper Pulp
Asphalt (SHRP)	Gums	Plastisols
Ceramic Slurries	Inks(screen printing)	Starches
Cosmetics	Organisols	Surface Coatings
Creams	Paints	Toothpaste
Dairy Products	Paper Coatings	Varnish

### HA/HB SERIES - HIGH VISCOSITY

Asphalt	Pastes
Caulking Compounds	Peanut Butter
Chocolate	Putty
Epoxies	Roofing Compounds
Gels	Sealants
Inks(ballpoint, offset, lithographic)	Sheet Molding Compound
Molasses	Tars



## Viscosity

Texture

Rheology

In-line

QC

Powder

Moisture

Gas



TEXTURE

RHEOLOGY

IN-LINE

VISCOSITY

POWDER

MOISTURE

GAS