

Dispermill Rotor-stator

Accessoires

The Rotor-stator is developed for emulsifying and homogenizing. Due to the powerful created turbulence, the dispersion process will take place in a little concentrated area where particles will decrease in size. By use of emulsifiers, the Rotor-stator can mix two products that are difficult to mix with usual dispersion methods.

The Rotor-stator shaft is interchangeable with your Discovery or X-proof shaft and comes with flange included. There are two versions available. The short version can be connected on the Discovery 100, and the X-Proof FR 075. The long version can be connected to the Discovery 200, 300 and the X-Proof FR 150.

Spare part	Order code
Rotor	21-1181
Stator	21-1182
Bushing	23-1183



FEATURES

- ✓ 40 mm Rotor-stator
- ✓ For emulsifying and homogenizing
- ✓ Made of stainless steel
- ✓ Interchangeable shaft



Dispermill Rotor-stator	Order code
Rotor-stator (Discovery 100, FR 075)	12-1008
Rotor-stator (Discovery 200, 300, FR 150)	12-1032